

implantswiss.com

NOVODENT SA
Y-Parc Technopole
Avenue des Sciences 11
CH-1400, Yverdon-les-Bains
Switzerland

 **implant**swiss
SMART
CHOICE

IMPLANT TECHNOLOGY

» Content

Implantswiss Dental System	2
Swiss Precision	4
Product Traceability	5
Glass Package	6
Platform Color Coding	7
» Implantswiss System Specifications	8
S.R.A. Surface	12
Hybrid Design	13
Conical Connection with Hexagonal Index	14
Platform Switching	14
Self Cutting Edge	15
Hygienic Neck Design	15
» Implantswiss Implants	16
Bone Level Implant	18
Tissue Level Implant	19
» Surgical Protocols	20
Implantswiss Surgical Drills	22
» Short Implant	28
Short Implant	29
Short Implant Surgical Protocol	31
Short Implant Surgical Drills	32
» Zygoma Implant	36
Zygoma Implant Surgical Protocol	39
Zygoma Implant Surgical Drills	39
» Implantswiss Prosthetic Components	40
» Bone Level System Prosthetic Components	42
Bone Level Healing Abutments	44
Couple System	44
Solid System	47
Octa System	48
O-ring System	49
Locator System	50
Multi-Unit System	51
Superb Overdenture System	52
Bar System	52
» Tissue Level System Prosthetic Components	54
Tissue Level Healing Abutments	56
Couple System	56
Solid System	58
Octa System	59
O-ring System	60
Locator System	61
Superb Overdenture System	62
» Zygoma Dental Implant System Prosthetic Components	63
» Drivers	64
» Screws	65
» Code Index	66

implantswiss

Switzerland

is at the heart of

science, mechanical design and craftsmanship



The well respected and acknowledged world authority...

Switzerland

has been one of the most important centers for craftsmanship that requires

high precision and expertise

throughout history.

The products with impeccable mechanical designs, can be used for a lifetime and add value to your life.

Implantswiss Dental System owes its **excellent** structure to the inspiration, experience and value from history.

Implantswiss is a strong and reliable brand in the dental industry with its innovative solutions, unique surface design and high quality.

Implantswiss is devoted to Research and Development through cooperations with universities and crucial investments in clinical research. Implantswiss follows the latest trends and initiates extensive and contemporary projects that may lead the future of implant dentistry.

Implantswiss Dental System


is the most advanced stage
in implant dentistry



SWISS PRECISION




Implantswiss products are manufactured with the state-of-the-arts technologies and high precision. Implantswiss never compromises on its quality and reliability; the success and high quality of Implantswiss lies in this perfectionist approach in each step of the manufacturing process.



PRODUCT TRACEABILITY

Each Implantswiss product can be easily traced from production to delivery.

Each product package has a barcode label with a lot number created by the ERP software system. All steps from production to stock, delivery and after sales can be monitored with the barcode system; therefore, the follow-up process can be conducted flawlessly.



GLASS PACKAGE

GLASS PACKAGE

The nature and human health have always been the primary concern for Implantswiss. Implantswiss, gives great importance to protection of nature and human health. Therefore, all Implantswiss Implants are packed in 100% natural and hygienic glass vials.

Since the glass structure doesn't contain any foreign material, glass packages are especially preferred for preserving the high value products.

The high barrier property, ease of sterilisation and unreactivity in chemical substances are among the many advantages of glass packages.

SMART PACKAGE


The implant is simply removed from the glass bottle and placed in the osteotomy using the fixture driver.

Each product package has a label containing all necessary information such as the type, the length and the color code of the product.

The color coding and associated diameters of the implants are classified as given.



PLATFORM COLOR CODING




TISSUE LEVEL

BONE LEVEL



» **IMPLANTSWISS
SYSTEM
SPECIFICATIONS**




» SRA* Surface


Following calcium phosphate blasting, the surface of the implants is etched using Double Acid Etching.

Implantswiss S.R.A. surface has a homogenous surface roughness that perfectly mimics the bone morphology.

* Sandblasted Rough Acid-etched Surface




3D-Area Parameters



Detail of the thread Bank and representation of the surface shape after removal of a 4th order polynomial

ISO 25178				
Height Parameters				
Sa	1,80	µm		
Sq	2,44	µm		
Feature Parameters				
S10z	15,1	µm	pruning = 5%	Ten point height
S5p	8,59	µm	pruning = 5%	Five point peak height
S5v	6,49	µm	pruning = 5%	Five point of height




Detail of the thread valley and representation of the surface shape after removal of a 4th order polynomial

ISO 25178				
Height Parameters				
Sa	1,25	µm		
Sq	1,60	µm		
Feature Parameters				
S10z	9,16	µm	pruning = 5%	Ten point height
S5p	3,39	µm	pruning = 5%	Five point peak height
S5v	5,76	µm	pruning = 5%	Five point of height

» Tapered form + Straight form = **Hybrid Design**

It is a unique implant design with cylindrical (straight) design on the neck area and conical (tapered) design through the apex of the implant.



The first 3 mm of Implantswiss implants running straight downward from the neck constitutes the straight form. The remaining portion constitutes the tapered form at a certain angle, depending on the length of the implant. The optimum combination of these two specifications creates the hybrid form. One advantage of the hybrid


design is to exert minimal stress on crestal cortical bone because of the straight portion at the neck, leading to an optimal load absorption rate. Due to its tapered form, 1/3 of the implant length is directly positioned in the osteotomy, which facilitate easy insertion of the implant in place.

» Conical connection with hexagonal index

creates an excellent implant-abutment connection.

Implantswiss implants create an excellent implant-abutment interface, which prevents the formation of micro-gap and micro-leakage thanks to the universal conical connection and hexagonal index.

In Implantswiss Dental Implant System, all abutments are compatible with all fixtures with any diameter (except for Q 3.3) making a user friendly implant system for dentists.



» Platform Switching

protects the alveolar bone around the neck area.

Micro-gap can be formed at the implant-abutment interface in Bone Level implants especially if they have a butt joint connection. Bacteria and by-products in these areas can cause serious problems such as bone resorption and periimplantitis.

Even though Implantswiss Implants have bacteria sealed conical connection, Switch Platform system further protects cervical bone by moving the implant-abutment interface away from alveolar bone. Implantswiss Implants with Switch Platform System are safer against bacteria and microleakage and have a longer life span.



» Self Cutting Edge

Implantswiss implants have a thread design which does not transfer stress to the bone during insertion. The cutting edges and self cutting feature provide ideal primary stability in any bone type without transferring stress thanks to its semi-aggressive thread design.

Unique thread design that does not transfer stress to the bone.




» Hygienic Neck Design

There are no micro threads on the neck area of Implantswiss implants therefore it is easier to maintain the hygiene when required.

Bacteria can easily colonize in the micro threads. In case of supra-alveolar positioning of the implant neck, the bone resorption risk is therefore minimized.

Easy to maintain the hygiene



» **IMPLANTSWISS
IMPLANTS**



BONE LEVEL IMPLANT



Bone Level implants are designed to provide excellent outcome in any bone quality and in any location.


Implantswiss bone level implants allow to design an ideal smile in the aesthetic zone.

	$\varnothing 3.3$	$\varnothing 3.7$	$\varnothing 4.3$	$\varnothing 4.8$	$\varnothing 5.5$
Length	$\varnothing 3.3$	$\varnothing 3.7$	$\varnothing 4.3$	$\varnothing 4.8$	$\varnothing 5.5$
8.0 mm	S-BMFSR3308	S-BFHR3708	S-BFHR4308	S-BFHR4808	S-BWPHR5508
10.0 mm	S-BMFSR3310	S-BFHR3710	S-BFHR4310	S-BFHR4810	S-BWPHR5510
12.0 mm	S-BMFSR3312	S-BFHR3712	S-BFHR4312	S-BFHR4812	S-BWPHR5512
14.0 mm	S-BMFSR3314	S-BFHR3714	S-BFHR4314	S-BFHR4814	-

TISSUE LEVEL IMPLANT


An excellent soft and hard tissue protection can be achieved with the smooth neck and location of implant-abutment interface of the Tissue Level implant.

The implant is designed to provide excellent outcomes in any case where the aesthetic outcome is not the primary concern.



» SURGICAL PROTOCOLS

BONE LEVEL - TISSUE LEVEL



Pilot drill

Extra sharp self stoppered drills can be used at low speed (300 rpm), and collect autogenous bone graft when used below 50 rpm.

Torque controlled Ratchet

Depth gauge

Point drill

Drill extension

Neck Drill

Hex Screw Driver adaptable with ratchet and hand piece

Parallel pins

Long and short implant drivers for Ø 3.3 Bone and Tissue Level systems adaptable with both ratchet and hand piece


IMPLANTSWISS FIXTURE DRILLS




SURGICAL DRILLS

Self stoppered surgical drills eliminate the risk of complications both during the preparation of the osteotomy and insertion of the implant. Implantswiss drills are recommended to be used at 200-300 rpm to prevent over heating and reduce mechanical trauma, which will shorten the osseo-integration time of the implant. The drills are designed to harvest autogenous bone when used below 50rpm without irrigation to be used as a graft material in regenerative procedures.

Implantswiss Surgical Kit contains a specific drill for each implant available in the system.




Implantswiss Surgical Kit includes various drills with different diameters and lengths. These drills are designed to simplify the surgical procedures resulting in shortened surgery time and decreased stress for the patients.




**ø3.3
Fixture**




D1 Bone



D2-3 Bone




D4 Bone




**ø3.7
Fixture**




D1 Bone



D2-3 Bone




D4 Bone




**ø4.3
Fixture**




D1 Bone



D2-3 Bone




D4 Bone




**ø4.8
Fixture**




D1 Bone



D2-3 Bone



D4 Bone



**ø5.5
Fixture**



D1 Bone

Point drill	2.2 drill	Pilot drill	3.3 Fixture drill	3.7 Fixture drill	4.3 Fixture drill	4.8 Fixture drill	5.5 Fixture drill	5.5 Neck drill	5.5 Fixture
-------------	-----------	-------------	-------------------	-------------------	-------------------	-------------------	-------------------	----------------	-------------

D2-3 Bone

Point drill	2.2 drill	Pilot drill	3.3 Fixture drill	3.7 Fixture drill	4.3 Fixture drill	4.8 Fixture drill	5.5 Fixture drill	5.5 Fixture
-------------	-----------	-------------	-------------------	-------------------	-------------------	-------------------	-------------------	-------------

D4 Bone

Point drill	2.2 drill	Pilot drill	3.3 Fixture drill	3.7 Fixture drill	4.3 Fixture drill	4.8 Fixture drill	5.5 Fixture
-------------	-----------	-------------	-------------------	-------------------	-------------------	-------------------	-------------

» **IMPLANTSWISS
SHORT IMPLANTS**




SHORT IMPLANT


In cases of insufficient bone volume, either above the mandibular canal or below the maxillary sinuses, short implants may be an alternative treatment to advanced surgical interventions aiming to regenerate hard tissues.

Due to aggressive thread design, Bone level short implants can be preferred mainly in maxilla, where the alveolar bone is relatively soft, whereas Tissue level short implant can be a good choice especially for the lower jaw.

BONE LEVEL




TISSUE LEVEL






SHORT IMPLANT SURGICAL PROTOCOL

The surgical kit is recommended to be used with both Bone and Tissue Level short implants.



SHORT IMPLANT DRILLS

Implantswiss Short Implant Surgical Kit contains a specific drill for each implant available in the system.

	Fixture Drills		
SURGICAL DRILLS	4.3	4.8	5.5
			

Implantswiss Short Implant Surgical Kit includes various drills with different diameters and lengths. These drills are designed to simplify the surgical procedures resulting in shortened surgery time and decreased stress for the patients.


AUXILIARY DRILLS	Point drill	Pilot drill
		

**ø4.3
Fixture**




BONE LEVEL

D4 / D3 Bone



TISSUE LEVEL

D2 / D1 Bone




Bone Level Short Implants are recommended to be used only in D3/D4 and Tissue Level Implants only in D1/D2 bone.

**ø4.8
Fixture**




BONE LEVEL

D4 / D3 Bone



TISSUE LEVEL

D2 / D1 Bone




**ø5.5
Fixture**




BONE LEVEL

D4 / D3 Bone



TISSUE LEVEL

D2 / D1 Bone



» **IMPLANTSWISS**
ZYGOMA IMPLANTS



ZYGOMA IMPLANT

Zygomatic implants are designed for cases where traditional implants can not be used and for patients who want to benefit from immediate loading protocols even in the presence of minimum residual maxillary bone height.

Implantswiss Zygomatic Implant offers three different product options.

ZYGOMA IMPLANT

ø4.3



35 mm	S-ZYG35
37.5 mm	S-ZYG37.5
40 mm	S-ZYG40
42.5 mm	S-ZYG42.5
45 mm	S-ZYG45
47.5 mm	S-ZYG47.5
50 mm	S-ZYG50
52.5 mm	S-ZYG52.5
55 mm	S-ZYG55
57.5 mm	S-ZYG57.5
60 mm	S-ZYG60

ATROPHIC MAXILLARY ZYGOMA IMPLANT

ø4.3



35 mm	S-ZAM35
37.5 mm	S-ZAM37.5
40 mm	S-ZAM40
42.5 mm	S-ZAM42.5
45 mm	S-ZAM45
47.5 mm	S-ZAM47.5
50 mm	S-ZAM50
52.5 mm	S-ZAM52.5
55 mm	S-ZAM55
57.5 mm	S-ZAM57.5
60 mm	S-ZAM60

ONCOLOGICAL ZYGOMA IMPLANT


ø4.3



25 mm	S-ONC25
27.5 mm	S-ONC27.5
30 mm	S-ONC30
32.5 mm	S-ONC32.5



ZYGOMA IMPLANT SURGICAL PROTOCOL



SURGICAL DRILLS

Point drill Pilot drill 3.5 Pilot drill Long Twist drill 2.5 Short Twist drill 2.5 Long Twist drill 2.9 Short Twist drill 2.9 Long Twist drill 3.5 Short Twist drill 3.5 Long


















» **IMPLANTSWISS
PROSTHETIC COMPONENTS**


















BONE LEVEL
















COUPLE SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	LAB ANALOG	ABUTMENT	SCREW
 <p>ø3.3 ø3.7 ø4.3</p>  <p>ø4.8 ø5.5</p>	 <p>Closing Screw</p>  <p>Healing Abutment</p>	<p>Pick-up</p>  <p>Hex N-Hex</p> <p>Transfer</p>  <p>Hex N-Hex</p>		<p>Couple</p>  <p>Hex N-Hex</p> <p>Multi</p>  <p>Hex N-Hex</p> <p>Angled</p>  <p>15° 25°</p>		
TOOL	 <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	 <p>Hex Screw Driver</p>			 <p>Hex Screw Driver</p>











SOLID SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	IMPRESSION CAP	LAB ANALOG	PLASTIC COPING
 <p>ø3.3 ø3.7 ø4.3</p>  <p>ø4.8 ø5.5</p>	 <p>Closing Screw</p>  <p>Healing Abutment</p>	<p>Solid</p>  <p>Solid</p> <p>Cap</p>  <p>Cap</p>			<p>Single</p>  <p>Single</p> <p>Bridge</p>  <p>Bridge</p>	
TOOL	 <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	 <p>Solid Driver</p>  <p>Hex Screw Driver</p>			 <p>Hex Screw Driver</p>

OCTA SYSTEM





FLOW	FIXTURE	HEALING	SUB OCTA ABUTMENT	IMPRESSION COPING	LAB ANALOG	LAB CYLINDER	SCREW
 <p>ø3.3 ø3.7 ø4.3</p>  <p>ø4.8 ø5.5</p>	 <p>Closing Screw</p>  <p>Healing Abutment</p>	 <p>Octa Healing Cap</p>  <p>Octa Abutment</p>	 <p>Pick-up</p>  <p>Transfer</p>		 <p>Octa</p> <p>N-Octa</p>		
TOOL	 <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	 <p>Octa Driver</p>				 <p>Hex Screw Driver</p>

O-RING SYSTEM

FLOW	FIXTURE	HEALING	O-RING ABUTMENT	O-RING LAB ANALOG	RETAINER
 <p>ø3.3 ø3.7 ø4.3</p>  <p>ø4.8 ø5.5</p>	 <p>Closing Screw</p>  <p>Healing Abutment</p>			 <p>O-Ring</p> <p>O-Ring Retainer</p>	
TOOL	 <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	 <p>O-Ring Driver</p>		




BONE LEVEL HEALING ABUTMENTS

BONE LEVEL HEALING ABUTMENTS

	 Ø3.7 <small>only for ø3.3 implant</small>	 Ø4.5	 Ø5.5	 Ø6.5
1 mm	-	S-BHA4510	S-BHA5510	S-BHA6510
2 mm	-	S-BHA4520	S-BHA5520	S-BHA6520
2.5 mm	S-BMHA3725	-	-	-
3 mm	-	S-BHA4530	S-BHA5530	S-BHA6530
3.5 mm	S-BMHA3735	-	-	-
4 mm	-	S-BHA4540	S-BHA5540	S-BHA6540
5 mm	S-BMHA3750	S-BHA4550	S-BHA5550	S-BHA6550
7 mm	S-BMHA3770	-	-	-

System » Bone Level
Material » Ti Grade 5 ELI
Driver » Hex Screw Driver
Torque Value » 15-20 Ncm





BONE LEVEL SHORT IMPLANT HEALING ABUTMENTS

	 Ø5.0	 Ø6.0	 Ø8.0
2 mm	S-BMHA5020	S-BMHA6020	S-BMHA8020
3 mm	S-BMHA5030	S-BMHA6030	S-BMHA8030

System » Bone Level
Material » Ti Grade 5 ELI
Driver » Hex Screw Driver
Torque Value » 15-20 Ncm




BONE LEVEL COUPLE SYSTEM

BONE LEVEL COUPLE ABUTMENTS

	 Ø3.7 <small>only for ø3.3 implant</small>	 Ø4.5	 Ø5.5	 Ø6.5
1 mm	S-BMCA37107	S-BCA45107	S-BCA55107	S-BCA65107
2 mm	S-BMCA37207	S-BCA45207	S-BCA55207	S-BCA65207
3 mm	S-BMCA37307	S-BCA45307	S-BCA55307	S-BCA65307
4 mm	S-BMCA37407	S-BCA45407	S-BCA55407	S-BCA65407
5 mm	S-BMCA37507	S-BCA45507	S-BCA55507	S-BCA65507





System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Hex Screw Driver

BONE LEVEL SHORT IMPLANT COUPLE ABUTMENTS

			
	ø5.0	ø6.0	ø8.0
2 mm	S-BMCA5020	S-BMCA6020	S-BMCA8020
3 mm	S-BMCA5030	S-BMCA6030	S-BMCA8030





System » Bone Level
 Material » Ti Grade 5 ELI
 Torque Value » 25 Ncm
 Driver » Hex Screw Driver

BONE LEVEL COUPLE ABUTMENTS NON-HEX

				
	ø3.7 <small>only for ø3.3 implant</small>	ø4.5	ø5.5	ø6.5
1 mm	S-BMCAN37107	S-BCAN45107	S-BCAN55107	S-BCAN65107
2 mm	S-BMCAN37207	S-BCAN45207	S-BCAN55207	S-BCAN65207
3 mm	S-BMCAN37307	S-BCAN45307	S-BCAN55307	S-BCAN65307
4 mm	S-BMCAN37407	S-BCAN45407	S-BCAN55407	S-BCAN65407
5 mm	S-BMCAN37507	S-BCAN45507	S-BCAN55507	S-BCAN65507





System » Bone Level
 Material » Ti Grade 5 ELI
 Torque Value » 25 Ncm
 Driver » Hex Screw Driver

BONE LEVEL ANGLED ABUTMENT

				
	ø3.7 <small>only for ø3.3 implant</small>	ø4.5 15° 25°	ø5.5 15° 25°	ø6.5 15° 25°
1 mm (15°)	S-BMAA37151	2 mm 15°	S-BAA45152	S-BAA55152
2 mm (15°)	S-BMAA37152	2 mm 25°	S-BAA45252	S-BAA55252
4 mm (15°)	S-BMAA37154	4 mm 15°	S-BAA45154	S-BAA55154
		4 mm 25°	S-BAA45254	S-BAA55254
			S-BAA65152	S-BAA65252
			S-BAA65154	S-BAA65254



System » Bone Level // Material » Ti Grade 5 ELI // Torque Value » 25 Ncm // Driver » Hex Screw Driver

BONE LEVEL ANGLED ABUTMENT NON-HEX

				
	ø3.7 <small>only for ø3.3 implant</small>	ø4.5 15° 25°	ø5.5 15° 25°	ø6.5 15° 25°
1 mm (15°)	S-BMAAN37151	2 mm 15°	S-BAAN45152	S-BAAN55152
2 mm (15°)	S-BMAAN37152	2 mm 25°	S-BAAN45252	S-BAAN55252
4 mm (15°)	S-BMAAN37154	4 mm 15°	S-BAAN45154	S-BAAN55154
		4 mm 25°	S-BAAN45254	S-BAAN55254
			S-BAAN65152	S-BAAN65252
			S-BAAN65154	S-BAAN65254




System » Bone Level // Material » Ti Grade 5 ELI // Torque Value » 25 Ncm // Driver » Hex Screw Driver

BONE LEVEL MULTI ABUTMENTS // HEX / NON-HEX

			
	ø3.7	ø4.5	ø5.5
HEX	S-BMMA37	S-BMA45	S-BMA55
NON-HEX	S-BMMAN37	S-BMAN45	S-BMAN55





System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Hex Screw Driver

TEMPORARY ABUTMENTS

			
	ø3.7 <small>only for ø3.3 implant</small>	ø4.5	ø5.5
	S-BMTA37	S-BTA45	S-BTA55



System » Bone Level
Material » PEEK
Torque Value » 20 Ncm
Driver » Hex Screw Driver

BONE LEVEL IMPRESSION COPING PICK UP

				
	ø3.7 <small>only for ø3.3 implant</small>	ø4.5	ø5.5	ø6.5
HEX	S-BMICP37L S-BMICP37S	S-BICP45L S-BICP45S	S-BICP55L S-BICP55S	S-BICP65L S-BICP65S
NON-HEX	S-BMICPN37L S-BMICPN37S	S-BICPN45L S-BICPN45S	S-BICPN55L S-BICPN55S	S-BICPN65L S-BICPN65S

System » Bone Level
Material » Ti Grade 5 ELI

BONE LEVEL SHORT IMPLANT IMPRESSION COPING PICK UP

		
	ø5.0	ø6.0
	S-BMSICP50	S-BMSICP60

System » Bone Level
Material » Ti Grade 5 ELI

BONE LEVEL IMPRESSION COPING TRANSFER

				
	ø3.7 <small>only for ø3.3 implant</small>	ø4.5	ø5.5	ø6.5
HEX	S-BMICT37L S-BMICT37S	S-BICT45L S-BICT45S	S-BICT55L S-BICT55S	S-BICT65L S-BICT65S
NON-HEX	S-BMICTN37L S-BMICTN37S	S-BICTN45L S-BICTN45S	S-BICTN55L S-BICTN55S	S-BICTN65L S-BICTN65S

System » Bone Level
Material » Ti Grade 5 ELI

BONE LEVEL SHORT IMPLANT IMPRESSION COPING TRANSFER



S-BMSICT50



S-BMSICT60

System » Bone Level
Material » Ti Grade 5 ELI

BONE LEVEL FIXTURE LAB ANALOG



S-BMFLA
only for ø3.3 implant



S-BFLA

System » Bone Level
Material » Ti Grade 5 ELI

BONE LEVEL SOLID SYSTEM

BONE LEVEL SOLID ABUTMENTS



ø3.7
only for ø3.3 implant



ø4.5



ø5.5



ø6.5

1 mm	S-BMSA3710	S-BSA4510	S-BSA5510	S-BSA6510
2 mm	S-BMSA3720	S-BSA4520	S-BSA5520	S-BSA6520
3 mm	S-BMSA3730	S-BSA4530	S-BSA5530	S-BSA6530
4 mm	S-BMSA3740	S-BSA4540	S-BSA5540	S-BSA6540
5 mm	S-BMSA3750	S-BSA4550	S-BSA5550	S-BSA6550

System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Hex Screw Driver

BONE LEVEL SOLID ABUTMENT CAP



ø3.7
only for ø3.3 implant

S-BMSACap37



ø5.1

S-BSACap45



ø5.8

S-BSACap55



ø6.8

S-BSACap65

System » Bone Level Solid System
Material » POM

BONE LEVEL SOLID ABUTMENT IMPRESSION CAP



ø3.7
only for ø3.3 implant

S-BMSAICap37



ø4.5

S-BSAICap45



ø5.5

S-BSAICap55



ø6.5

S-BSAICap65

System » Bone Level Solid System
Material » Plastic

BONE LEVEL SOLID LAB ANALOG



ø3.7

only for ø3.3 implant

S-BMSALA37



ø4.5

S-BSALA45



ø5.5

S-BSALA55



ø6.5

S-BSALA65

System » Bone Level Solid System

Material » Ti Grade 5 ELI

BONE LEVEL OCTA SYSTEM

BONE LEVEL OCTA ABUTMENTS



ø3.7

only for ø3.3 implant

1 mm

S-BMOA10



ø4.8

S-BOA10

2 mm

S-BMOA20

S-BOA20

3 mm

S-BMOA30

S-BOA30

4 mm

S-BMOA40

S-BOA40

5 mm

S-BMOA50

S-BOA50

System » Bone Level

Material » Ti Grade 5 ELI

Torque Value » : Abutment Screw
25 Ncm Oclusal Screw 20 Ncm

Driver » Octa Driver

BONE LEVEL OCTA ABUTMENT HEALING CAP



ø4.8

4 mm

S-OAHCap

6 mm

S-OAHCapL

System » Octa System

Material » Ti Grade 5 ELI

BONE LEVEL OCTA SYSTEM PLASTIC CYLINDER



OCTA

S-OAPICy

NON-OCTA

S-OAPICyN

System » Octa System

Material » POM

BONE LEVEL OCTA SYSTEM TITANIUM CYLINDER

HEX

S-OATICy

NON-HEX

S-OATICyN

System » Octa System

Material » Ti Grade 5 ELI

BONE LEVEL OCTA SYSTEM IMPRESSION COPING PICK UP



ø5.5

HEX S-ICOAP
NON-HEX S-ICOAPN

BONE LEVEL OCTA SYSTEM IMPRESSION COPING TRANSFER



ø5.0

HEX S-ICOAT
NON-HEX S-ICOATN

System » Octa System
Material » Ti Grade 5 ELI

BONE LEVEL OCTA ABUTMENT LAB ANALOG



ø4.8

S-OALA

System » Octa System
Material » Ti Grade 5 ELI

BONE LEVEL O-RING SYSTEM

BONE LEVEL O-RING ABUTMENTS



ø3.42



ø4.5



ø2.9

only for ø3.3 implant



ø3.7

only for ø3.3 implant

	ø3.42	ø4.5	ø2.9 <small>only for ø3.3 implant</small>	ø3.7 <small>only for ø3.3 implant</small>
0.5 mm	S-BORA05	-	S-BMORA05	-
1 mm	S-BORA10	-	-	-
2 mm	-	S-BORA20	-	S-BMORA20
3 mm	-	S-BORA30	-	-
4 mm	-	S-BORA40	-	S-BMORA40
6 mm	-	S-BORA60	-	-

System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » O-Ring Driver

BONE LEVEL RETAINER



S-ORR



S-ORRc



S-ORRo

TITANIUM CAP

GREEN RETAINER
Retention 350 Gr



PINK RETAINER
Retention 900 Gr



YELLOW RETAINER
Retention 500 Gr

BONE LEVEL O-RING LAB ANALOG



S-ORLA

System » Bone Level O-Ring System
Material » Ti Grade 5 ELI

BONE LEVEL LOCATOR SYSTEM

BONE LEVEL MINI LOCATOR ABUTMENTS

ø3.3
only for ø3.3 implant

1.5 mm	S-BMLA15
2.5 mm	S-BMLA25
3.5 mm	S-BMLA35
4.5 mm	S-BMLA45
5.5 mm	S-BMLA55
6.5 mm	S-BMLA65

BONE LEVEL LOCATOR ABUTMENTS

ø3.86

1.5 mm	S-BLA15
2.5 mm	S-BLA25
3.5 mm	S-BLA35
4.5 mm	S-BLA45
5.5 mm	S-BLA55
6.5 mm	S-BLA65



System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Locator Driver

BONE LEVEL LOCATOR MATRIX



S-LMht



VIOLET CAP
Rigid Retention 2.7 Kg



YELLOW CAP
Extra Soft Retention 0.6 Kg



CLEAR CAP
Standard Retention 1.8 Kg



BLACK CAP
Processing



PINK CAP
Soft Retention 1.2 Kg

BONE LEVEL LOCATOR LAB ANALOG



S-BLLA Ø 4mm

BONE LEVEL MULTI UNIT SYSTEM

MULTI UNIT ABUTMENTS



∅4.8 (17°)

2.5 mm S-BMUA48172

3.5 mm S-BMUA48173



∅4.8 (35°)

S-BMUA48352

S-BMUA48353

System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » : Abutment Screw
25 Ncm Occlusal Screw 20 Ncm
Driver » 1,2 Hex Driver

MULTI UNIT ABUTMENT HEALING CAP



∅4.8

S-BMUAHCap

System » Multi Unit System
Material » Ti Grade 5 ELI

MULTI UNIT IMPRESSION COPING TRANSFER



∅5.0

NON-HEX S-ICBMUATN

System » Multi Unit System
Material » Ti Grade 5 ELI

MULTI UNIT ABUTMENT PLASTIC CYLINDER



S-MUAPICY

MULTI UNIT LAB ANALOG



S-MUALA
∅ 4.8×14mm

System » Bone Level
Material » Ti Grade 5 ELI

BONE LEVEL SUPERB OVERDENTURE SYSTEM

SUPERB OVERDENTURE ABUTMENT

	ø2.9	ø3.5	ø3.42	ø4.5
0.5 mm	S-BMBA05	-	S-BEBA05	-
1 mm	-	-	S-BEBA10	-
2 mm	-	S-BMEBA20	-	S-BEBA20
3 mm	-	-	-	S-BEBA30
4 mm	-	S-BMEBA40	-	S-BEBA40
5 mm	-	-	-	S-BEBA50
6 mm	-	-	-	S-BEBA60



System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Superb Overdenture Driver

SUPERB OVERDENTURE CAP



S-EBOCap

Ø 5.5mmx3.5mm

SUPERB OVERDENTURE LAB ANALOG



S-EBLA

Ø 4.1mmx16,6 mm

BAR SYSTEM

BAR BURNOUT

S-IB - burnout



BAR POSITIONING

S-IBP - burnout



BAR HOUSING

S-IBH - burnout



BAR CLIP YELLOW

S-IBClip - y

Medium Retention 1.8 Kg



BAR CLIP PINK

S-IBClip - p

Soft Retention 1.2 Kg



BAR CLIP RED

S-IBClip - r

Extra Soft Retention 0.8 Kg




TOOL FOR INSERTING CLIP



KEY FOR PARALLELOMETER



















SUPERB OVERDENTURE SYSTEM


















TISSUE LEVEL

COUPLE SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	LAB ANALOG	ABUTMENT	SCREW
 <p>ø3.7 ø4.3-L06 ø4.3</p>  <p>ø4.8-L06 ø4.8 ø5.5</p>	 <p>Cover Screw</p>  <p>Healing Abutment</p>	<p>Pick-up</p>  <p>Hex N-Hex</p> <p>Transfer</p>  <p>Hex N-Hex</p>		<p>Couple</p>  <p>Hex N-Hex</p> <p>Multi</p>  <p>Hex N-Hex</p> <p>Angled</p>  <p>15° 25°</p> <p>Synocta</p>  <p>Hex N-Hex</p>		
<p>TOOL</p>  <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	 <p>Hex Screw Driver</p>			 <p>Hex Screw Driver</p>	

SOLID SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	IMPRESSION CAP	LAB ANALOG	PLASTIC COPING
 <p>ø3.7 ø4.3-L06 ø4.3</p>  <p>ø4.8-L06 ø4.8 ø5.5</p>	 <p>Closing Screw</p>  <p>Healing Abutment</p>	  <p>Cap</p>		<p>Shoulder Analog & Pin</p>  <p>Lab Analog</p> 	<p>Crown</p>  <p>Bridge</p> 	
<p>TOOL</p>  <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	  <p>Solid Driver Hex Screw Driver</p>			 <p>Hex Screw Driver</p>	

OCTA SYSTEM

FLOW	FIXTURE	HEALING	OCTA ABUTMENT	IMPRESSION COPING	LAB ANALOG	LAB CYLINDER	SCREW
TOOL	 Ratchet Driver	 Hex Screw Driver	 Octa Driver				 Hex Screw Driver

O-RING SYSTEM

FLOW	FIXTURE	HEALING	O-RING ABUTMENT	O-RING LAB ANALOG	RETAINER
TOOL	 Ratchet Driver	 Hex Screw Driver	 O-Ring Driver		

TISSUE LEVEL HEALING ABUTMENTS

TISSUE LEVEL HEALING ABUTMENTS



ø3.7

ø6.4

3 mm	S-THA30	S-THAW30
4 mm	S-THA40	S-THAW40
5.5 mm	S-THA55	S-THAW55

System » Tissue Level
Material » Ti Grade 4
Torque Value » 20 Ncm
Driver » Hex Screw Driver

TISSUE LEVEL COUPLE SYSTEM

TISSUE LEVEL COUPLE ABUTMENTS

TISSUE LEVEL COUPLE ABUTMENT NON-OCTA



ø5.2

ø6.4

1 mm	S-TCA10	S-TCAW10
2 mm	S-TCA20	S-TCAW20
3 mm	S-TCA30	S-TCAW30
4 mm	S-TCA40	S-TCAW40

ø5.2

ø6.4

1 mm	S-TCAN10	S-TCAWN10
2 mm	S-TCAN20	S-TCAWN20
3 mm	S-TCAN30	S-TCAWN30
4 mm	S-TCAN40	S-TCAWN40

System » Tissue Level // **Material** » Ti Grade 4 // **Torque Value** » 25 Ncm // **Driver** » Hex Screw Driver

TISSUE LEVEL ANGLED ABUTMENT



ø3.7

15°	S-TAA15
25°	S-TAA25

System » Tissue Level
Material » Ti Grade 4
Torque Value » 25 Ncm
Driver » Hex Screw Driver

TISSUE LEVEL MULTI ABUTMENT OCTA / NON-OCTA



ø5.2

OCTA	S-TMA
N-OCTA	S-TMAN

System » Tissue Level
Material » Ti Grade 4
Torque Value » 25 Ncm
Driver » Hex Screw Driver

TISSUE LEVEL SYNOCTA ABUTMENT



ø3.5

6 mm

S-TSOA

TISSUE LEVEL SYNOCTA ABUTMENT NON-OCTA



ø3.5

6 mm

S-TSOAN

System » Tissue Level // Material » Ti Grade 4 // Torque Value » 25 Ncm // Driver » Hex Screw Driver

TISSUE LEVEL IMPRESSION COPING PICK UP OCTA / NON-OCTA



ø5.5

OCTA

S-TICP

N-OCTA

S-TICPN



ø6.4

S-TICPW

S-TICPWN

System » Tissue Level
Material » Ti Grade 5 ELI

TISSUE LEVEL IMPRESSION COPING TRANSFER OCTA / NON-OCTA



ø5.5

OCTA

S-TICT

N-OCTA

S-TICTN

System » Tissue Level
Material » Ti Grade 5 ELI

TISSUE LEVEL FIXTURE LAB ANALOG



ø4.8

S-TFLA



ø6.0

S-TFLLA

System » Tissue Level
Material » Ti Grade 5 ELI

TISSUE LEVEL SOLID SYSTEM

TISSUE LEVEL SOLID ABUTMENTS



ø3.5 ø4.5

4 mm	S-TSA40	S-TSAL40
5 mm	-	S-TSAL50
5.5 mm	S-TSA55	-
7 mm	S-TSA70	S-TSAL70

System » Tissue Level
Material » Ti Grade 4
Torque Value » 25 Ncm
Driver » Solid Driver

TISSUE LEVEL SOLID ABUTMENT POSITIONING CYLINDER



S-TSAPoCy

System » Tissue Level Solid System
Material » Plastic

TISSUE LEVEL SOLID PLASTIC COPING



ø5.3

CROWN	S-TSAPICc
BRIDGE	S-TSAPICb

System » Tissue Level Solid System
Material » POM

TISSUE LEVEL SOLID ABUTMENT IMPRESSION CAP



S-TSAICap

System » Tissue Level Solid System
Material » Plastic

TISSUE LEVEL SOLID ABUTMENT LAB ANALOG



ø4.8



ø6.0

4 mm	S-TSALA40	S-TSALLA40
5.5 mm	S-TSALA55	S-TSALLA55
7 mm	S-TSALA70	S-TSALLA70

System » Tissue Level Solid System
Material » Ti Grade 5 ELI

TISSUE LEVEL OCTA SYSTEM

TISSUE LEVEL OCTA ABUTMENT



ø3.5

S-TOA

System » Tissue Level
Material » Ti Grade 4
Torque Value » : Abutment Screw
25 Ncm Oclusal Screw 20 Ncm
Driver » Octa Driver

TISSUE LEVEL OCTA ABUTMENT HEALING CAP



ø4.8

4 mm	S-OAHCap
6 mm	S-OAHCapL

System » Octa System
Material » Ti Grade 5 ELI

TISSUE LEVEL OCTA ABUTMENT PLASTIC CYLINDER OCTA / NON-OCTA



OCTA	S-OAPICy
N-OCTA	S-OAPICyN

System » Octa System
Material » POM

TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING PICK UP OCTA / NON-OCTA



ø5.5

OCTA S-ICOAP
N-OCTA S-ICOAPN

System » Octa System
Material » Ti Grade 5 ELI

TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING TRANSFER OCTA / NON-OCTA



ø5.0

OCTA S-ICOAT
N-OCTA S-ICOATN

System » Octa System
Material » Ti Grade 5 ELI

TISSUE LEVEL OCTA LAB ANALOG



ø4.8

S-OALA

System » Octa System
Material » Ti Grade 5 ELI

TISSUE LEVEL O-RING SYSTEM

TISSUE LEVEL RETAINER



S-ORR



S-ORRc



TITANIUM CAP



PINK RETAINER / Retention 900 Gr



GREEN RETAINER / Retention 350 Gr



YELLOW RETAINER / Retention 500 Gr

S-ORRo

TISSUE LEVEL O-RING ABUTMENTS



Ø3.5

0.5 mm	S-TORA05
2 mm	S-TORA20
4 mm	S-TORA40

System » Tissue Level
Material » Ti Grade 4
Torque Value » 25 Ncm
Driver » O-Ring Driver

TISSUE LEVEL O-RING LAB ANALOG



Ø 4.1

S-ORLA

System » O-Ring System
Material » Ti Grade 5 ELI

TISSUE LEVEL LOCATOR SYSTEM

TISSUE LEVEL LOCATOR ABUTMENTS

Ø4.9

1.5 mm	S-TLA15
2.5 mm	S-TLA25
3.5 mm	S-TLA35
4.5 mm	S-TLA45
5.5 mm	S-TLA55
6.5 mm	S-TLA65



System » Tissue Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Locator Driver

TISSUE LEVEL LOCATOR MATRIX



S-LMht



VIOLET CAP
Rigid Retention 2.7 Kg



YELLOW CAP
Extra Soft Retention 0.6 Kg



CLEAR CAP
Standard Retention 1.8 Kg



BLACK CAP
Processing



PINK CAP
Soft Retention 1.2 Kg

TISSUE LEVEL LOCATOR LAB ANALOG



Ø 4.90

S-TLLA

TISSUE LEVEL SUPERB OVERDENTURE SYSTEM

TISSUE LEVEL SUPERB OVERDENTURE ABUTMENT

Ø3.5

0.5 mm S-TEBA05
2 mm S-TEBA20
4 mm S-TEBA40

System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Superb Overdenture Driver

TISSUE LEVEL SUPERB OVERDENTURE CAP



S-EBOCap

TISSUE LEVEL SUPERB OVERDENTURE LAB ANALOG



16,6 mm

Ø 4.1

S-EBLA

ZYGOMA DENTAL IMPLANT SYSTEM

ABUTMENTS

ABUTMENTS



3 mm S-ZCA4330
5 mm S-ZCA4350

ANGLED ABUTMENTS



3 mm S-ZAA43203
5 mm S-ZAA43205

System » Zygoma
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Hex Screw Driver

MULTI UNIT



17° S-ZMUA43172

System » Zygoma
Material » Ti Grade 5 ELI
Driver » Hex Screw Driver

LAB ANALOG



S-ZLA

System » Zygoma // **Material** » Ti Grade 5 ELI

CAPS

HEALING CAP



3 mm S-ZHA4330
5 mm S-ZHA4350

System » Zygoma // **Torque Value** » 25 Ncm

COPING

IMPRESSION COPING PICK UP



S-ZICP

System » Zygoma
Material » Ti Grade 5 ELI

PLASTIC CYLINDER



S-ZPICYN

System » Zygoma
Material » POM

DRIVERS

Implant Drivers (For ratchet and contra angle)

Bone Level Mini Fixture Driver Machine Short	S-BMFDMS
Bone Level Mini Fixture Driver Machine Long	S-BMFDML
Bone Level Mini Fixture Driver Ratchet Short	S-BMFDRS
Bone Level Mini Fixture Driver Ratchet Long	S-BMFDRL
Bone Level Fixture Driver Machine Short	S-BFDMS
Bone Level Fixture Driver Machine Long	S-BFDML
Bone Level Fixture Driver Ratchet Short	S-BFDRS
Bone Level Fixture Driver Ratchet Long	S-BFDRL
Tissue Level Fixture Driver Ratchet Short	S-TFDRS
Tissue Level Fixture Driver Ratchet Long	S-TFDRL
Tissue Level Fixture Driver Machine Short	S-TFDMS
Tissue Level Fixture Driver Machine Long	S-TFDML
Hex driver (extra short, short, long, extra long)	S-HDRS / S-HDRL
Hex Driver Contra Angle (extra short, short, long, extra long)	S-HDRS / S-HDRL
Superb Overdenture Driver	S-EBRD
O-ring Driver	S-ORADS/L
Octa Driver	S-OADS/L
Solid Driver	S-SADL
ORRO Insert Adaptor	S-OIA
Bar Matrix Insert Adaptor	S-BMIA

SCREWS

CLOSING SCREW

BONE LEVEL

Bone Level Mini Closing Screw	ø 2.8	4.5 mm	S-BMCS
Bone Level Closing Screw	ø 3.36	6.5 mm	S-BCS

TISSUE LEVEL

Tissue Level Closing-Screw Small	ø 3.5	5.35 mm	S-TCSS
Tissue Level Closing-Screw Large	ø 5.2	6.4mm	S-TCSL
Tissue Level Closing-Screw Wide	ø 6.4	6.4 mm	S-TCSW

ABUTMENT SCREWS

BONE LEVEL

Bone Level Mini Abutment Screw	ø 2.1	7.64 mm	S-BMAS
Bone Level Abutment Screw	ø 2.3	9.32 mm	S-BAS
Bone Level Multi-Unit Abutment Screw	ø 1.96	9 mm	S-BMUAS

TISSUE LEVEL

Tissue Level Abutment Screw	ø 2.5	7.1	S-TAS
-----------------------------	-------	-----	-------

COMMON SCREWS

Octa Abutment Screw	ø 2.5	4.5 mm	S-OAS
Multi-Unit Abutment Healing Screw	ø 2.2	4 mm	S-MUAHS

INDEX CODE

PRODUCT	PLATFORM	HEIGHT	CODE
BONE LEVEL FIXTURES			
BONE LEVEL MINI FIXTURE			
Bone Level Mini Fixture Straight SRA	ø 3.3	8 mm	S-BMFSR3308
Bone Level Mini Fixture Straight SRA	ø 3.3	10 mm	S-BMFSR3310
Bone Level Mini Fixture Straight SRA	ø 3.3	12 mm	S-BMFSR3312
Bone Level Mini Fixture Straight SRA	ø 3.3	14 mm	S-BMFSR3314
BONE LEVEL SHORT FIXTURE			
Bone Level Short Fixture Straight SRA	ø 4.3	5 mm	S-BMFSR4305
Bone Level Short Fixture Straight SRA	ø 4.3	6 mm	S-BMFSR4306
Bone Level Short Fixture Straight SRA	ø 4.8	5 mm	S-BMFSR4805
Bone Level Short Fixture Straight SRA	ø 4.8	6 mm	S-BMFSR4806
Bone Level Short Fixture Straight SRA	ø 5.5	5 mm	S-BMFSR5505
Bone Level Short Fixture Straight SRA	ø 5.5	6 mm	S-BMFSR5506
BONE LEVEL FIXTURE			
Bone Level Wide Fixture Hybrid SRA	ø 3.7	8 mm	S-BFHR3708
Bone Level Wide Fixture Hybrid SRA	ø 3.7	10 mm	S-BFHR3710
Bone Level Wide Fixture Hybrid SRA	ø 3.7	12 mm	S-BFHR3712
Bone Level Wide Fixture Hybrid SRA	ø 3.7	14 mm	S-BFHR3714
Bone Level Wide Fixture Hybrid SRA	ø 4.3	8 mm	S-BFHR4308
Bone Level Wide Fixture Hybrid SRA	ø 4.3	10 mm	S-BFHR4310
Bone Level Wide Fixture Hybrid SRA	ø 4.3	12 mm	S-BFHR4312
Bone Level Wide Fixture Hybrid SRA	ø 4.3	14 mm	S-BFHR4314
Bone Level Wide Fixture Hybrid SRA	ø 4.8	8 mm	S-BFHR4808
Bone Level Wide Fixture Hybrid SRA	ø 4.8	10 mm	S-BFHR4810
Bone Level Wide Fixture Hybrid SRA	ø 4.8	12 mm	S-BFHR4812
Bone Level Wide Fixture Hybrid SRA	ø 4.8	14 mm	S-BFHR4814
BONE LEVEL WIDE FIXTURE			
Bone Level Wide Fixture Hybrid SRA	ø 5.5	8 mm	S-BFHR5508
Bone Level Wide Fixture Hybrid SRA	ø 5.5	10 mm	S-BFHR5510
Bone Level Wide Fixture Hybrid SRA	ø 5.5	12 mm	S-BFHR5512
HEALING ABUTMENTS			
HEALING ABUTMENTS			
Healing Abutment	ø 3.5	2.5 mm	S-BMHA3725
Healing Abutment	ø 3.5	3.5 mm	S-BMHA3735
Healing Abutment	ø 3.5	6 mm	S-BMHA3750
Healing Abutment	ø 3.5	7 mm	S-BMHA3770
Healing Abutment	ø 4.5	1 mm	S-BHA4510
Healing Abutment	ø 4.5	2 mm	S-BHA4520

PRODUCT	PLATFORM	HEIGHT	CODE
Healing Abutment	∅ 4.5	3 mm	S-BHA4530
Healing Abutment	∅ 4.5	4 mm	S-BHA4540
Healing Abutment	∅ 4.5	5 mm	S-BHA4550
Healing Abutment	∅ 5.5	1 mm	S-BHA5510
Healing Abutment	∅ 5.5	2 mm	S-BHA5520
Healing Abutment	∅ 5.5	3 mm	S-BHA5530
Healing Abutment	∅ 5.5	4 mm	S-BHA5540
Healing Abutment	∅ 5.5	5 mm	S-BHA5550
Healing Abutment	∅ 6.5	1 mm	S-BHA6510
Healing Abutment	∅ 6.5	2 mm	S-BHA6520
Healing Abutment	∅ 6.5	3 mm	S-BHA6530
Healing Abutment	∅ 6.5	4 mm	S-BHA6540
Healing Abutment	∅ 6.5	5 mm	S-BHA6550
Mini Healing Abutment for Short Implant	∅ 5.0	2 mm	S-BMHA5020
Mini Healing Abutment for Short Implant	∅ 5.0	3 mm	S-BMHA5030
Mini Healing Abutment for Short Implant	∅ 6.0	2 mm	S-BMHA6020
Mini Healing Abutment for Short Implant	∅ 6.0	3 mm	S-BMHA6030
Mini Healing Abutment for Short Implant	∅ 8.0	2 mm	S-BMHA8020
Mini Healing Abutment for Short Implant	∅ 8.0	3 mm	S-BMHA8030

BONE LEVEL COUPLE SYSTEM

COUPLE ABUTMENTS

Mini Couple Abutment H7	∅ 3.5	1 mm	S-BMCA37107
Mini Couple Abutment H7	∅ 3.5	2 mm	S-BMCA37207
Mini Couple Abutment H7	∅ 3.5	3 mm	S-BMCA37307
Mini Couple Abutment H7	∅ 3.5	4 mm	S-BMCA37407
Mini Couple Abutment H7	∅ 3.5	5 mm	S-BMCA37507
Couple Abutment H7	∅ 4.5	1 mm	S-BCA45107
Couple Abutment H7	∅ 4.5	2 mm	S-BCA45207
Couple Abutment H7	∅ 4.5	3 mm	S-BCA45307
Couple Abutment H7	∅ 4.5	4 mm	S-BCA45407
Couple Abutment H7	∅ 4.5	5 mm	S-BCA45507
Couple Abutment H7	∅ 5.5	1 mm	S-BCA55107
Couple Abutment H7	∅ 5.5	2 mm	S-BCA55207
Couple Abutment H7	∅ 5.5	3 mm	S-BCA55307
Couple Abutment H7	∅ 5.5	4 mm	S-BCA55407
Couple Abutment H7	∅ 5.5	5 mm	S-BCA55507
Couple Abutment H7	∅ 6.5	1 mm	S-BCA65107
Couple Abutment H7	∅ 6.5	2 mm	S-BCA65207
Couple Abutment H7	∅ 6.5	3 mm	S-BCA65307
Couple Abutment H7	∅ 6.5	4 mm	S-BCA65407
Couple Abutment H7	∅ 6.5	5 mm	S-BCA65507

PRODUCT	PLATFORM	HEIGHT	CODE
Mini Couple Abutment for Short Implant	ø 5.0	2 mm	S-BMCA5020
Mini Couple Abutment for Short Implant	ø 5.0	3 mm	S-BMCA5030
Mini Couple Abutment for Short Implant	ø 6.0	2 mm	S-BMCA6020
Mini Couple Abutment for Short Implant	ø 6.0	3 mm	S-BMCA6030
Mini Couple Abutment for Short Implant	ø 8.0	2 mm	S-BMCA8020
Mini Couple Abutment for Short Implant	ø 8.0	3 mm	S-BMCA8030
COUPLE ABUTMENTS NON-HEX			
Mini Couple Abutment H7 Non-Hex	ø 3.5	1 mm	S-BMCAN3710
Mini Couple Abutment H7 Non-Hex	ø 3.5	2 mm	S-BMCAN3720
Mini Couple Abutment H7 Non-Hex	ø 3.5	3 mm	S-BMCAN3730
Mini Couple Abutment H7 Non-Hex	ø 3.5	4 mm	S-BMCAN3740
Mini Couple Abutment H7 Non-Hex	ø 3.5	5 mm	S-BMCAN3750
Couple Abutment H7 Non-Hex	ø 4.5	1 mm	S-BCAN4510
Couple Abutment H7 Non-Hex	ø 4.5	2 mm	S-BCAN4520
Couple Abutment H7 Non-Hex	ø 4.5	3 mm	S-BCAN4530
Couple Abutment H7 Non-Hex	ø 4.5	4 mm	S-BCAN4540
Couple Abutment H7 Non-Hex	ø 4.5	5 mm	S-BCAN4550
Couple Abutment H7 Non-Hex	ø 5.5	1 mm	S-BCAN5510
Couple Abutment H7 Non-Hex	ø 5.5	2 mm	S-BCAN5520
Couple Abutment H7 Non-Hex	ø 5.5	3 mm	S-BCAN5530
Couple Abutment H7 Non-Hex	ø 5.5	4 mm	S-BCAN5540
Couple Abutment H7 Non-Hex	ø 5.5	5 mm	S-BCAN5550
Couple Abutment H7 Non-Hex	ø 6.5	1 mm	S-BCAN6510
Couple Abutment H7 Non-Hex	ø 6.5	2 mm	S-BCAN6520
Couple Abutment H7 Non-Hex	ø 6.5	3 mm	S-BCAN6530
Couple Abutment H7 Non-Hex	ø 6.5	4 mm	S-BCAN6540
Couple Abutment H7 Non-Hex	ø 6.5	5 mm	S-BCAN6550
ANGLED ABUTMENTS			
Mini Angled Abutment	ø 3.5	1 mm 15°	S-BMAA37151
Mini Angled Abutment	ø 3.5	2 mm 15°	S-BMAA37152
Mini Angled Abutment	ø 3.5	4 mm 15°	S-BMAA37154
Angled Abutment	ø 4.5	2 mm 15°	S-BAA45152
Angled Abutment	ø 4.5	2 mm 25°	S-BAA45252
Angled Abutment	ø 4.5	4 mm 15°	S-BAA45154
Angled Abutment	ø 4.5	4 mm 25°	S-BAA45254
Angled Abutment	ø 5.5	2 mm 15°	S-BAA55152
Angled Abutment	ø 5.5	2 mm 25°	S-BAA55252
Angled Abutment	ø 5.5	4 mm 15°	S-BAA55154
Angled Abutment	ø 5.5	4 mm 25°	S-BAA55254
Angled Abutment	ø 6.5	2 mm 15°	S-BAA65152
Angled Abutment	ø 6.5	2 mm 25°	S-BAA65252
Angled Abutment	ø 6.5	4 mm 15°	S-BAA65154
Angled Abutment	ø 6.5	4 mm 25°	S-BAA65254

PRODUCT	PLATFORM	HEIGHT	CODE
ANGLED ABUTMENTS NON-HEX			
Mini Angled Abutment Non-Hex	ø 3.5	1 mm 15°	S-BMAAN37151
Mini Angled Abutment Non-Hex	ø 3.5	2 mm 15°	S-BMAAN37152
Mini Angled Abutment Non-Hex	ø 3.5	4 mm 15°	S-BMAAN37154
Angled Abutment Non-Hex	ø 4.5	2 mm 15°	S-BAAN45152
Angled Abutment Non-Hex	ø 4.5	2 mm 25°	S-BAAN45252
Angled Abutment Non-Hex	ø 4.5	4 mm 15°	S-BAAN45154
Angled Abutment Non-Hex	ø 4.5	4 mm 25°	S-BAAN45254
Angled Abutment Non-Hex	ø 5.5	2 mm 15°	S-BAAN55152
Angled Abutment Non-Hex	ø 5.5	2 mm 25°	S-BAAN55252
Angled Abutment Non-Hex	ø 5.5	4 mm 15°	S-BAAN55154
Angled Abutment Non-Hex	ø 5.5	4 mm 25°	S-BAAN55254
Angled Abutment Non-Hex	ø 6.5	2 mm 15°	S-BAAN65152
Angled Abutment Non-Hex	ø 6.5	2 mm 25°	S-BAAN65252
Angled Abutment Non-Hex	ø 6.5	4 mm 15°	S-BAAN65154
Angled Abutment Non-Hex	ø 6.5	4 mm 25°	S-BAAN65254
MULTI ABUTMENTS			
Mini Multi Abutment	ø 3.5		S-BMMA37
Multi Abutment	ø 4.5		S-BMA45
Multi Abutment	ø 5.5		S-BMA55
MULTI ABUTMENTS NON-HEX			
Mini Multi Abutment Non-Hex	ø 3.5		S-BMMAN37
Multi Abutment Non-Hex	ø 4.5		S-BMAN45
Multi Abutment Non-Hex	ø 5.5		S-BMAN55
TEMPORARY ABUTMENTS			
Mini Temporary Abutment	ø 3.7		S-BMTA37
Temporary Abutment	ø 4.5		S-BTA45
Temporary Abutment	ø 5.5		S-BTA55
IMPRESSION COPING PICK UP			
Mini Impression Coping Pick up	ø 3.5		S-BMICP37L
Mini Impression Coping Pick up	ø 3.5		S-BMICP37S
Impression Coping Pick up	ø 4.5		S-BICP45L
Impression Coping Pick up	ø 4.5		S-BICP45S
Impression Coping Pick up	ø 5.5		S-BICP55L
Impression Coping Pick up	ø 5.5		S-BICP55S
Impression Coping Pick up	ø 6.5		S-BICP65L
Impression Coping Pick up	ø 6.5		S-BICP65S
IMPRESSION COPING PICK UP NON-HEX			
Mini Impression Coping Pick up Non-Hex	ø 3.5		S-BMICPN37L
Mini Impression Coping Pick up Non-Hex	ø 3.5		S-BMICPN37S

PRODUCT	PLATFORM	HEIGHT	CODE
Impression Coping Pick up Non-Hex	ø 4.5		S-BICPN45L
Impression Coping Pick up Non-Hex	ø 4.5		S-BICPN45S
Impression Coping Pick up Non-Hex	ø 5.5		S-BICPN55L
Impression Coping Pick up Non-Hex	ø 5.5		S-BICPN55S
Impression Coping Pick up Non-Hex	ø 6.5		S-BICPN65L
Impression Coping Pick up Non-Hex	ø 6.5		S-BICPN65S
IMPRESSION COPING TRANSFER			
Mini Impression Coping Transfer	ø 3.5		S-BMICT37L
Mini Impression Coping Transfer	ø 3.5		S-BMICT37S
Impression Coping Transfer	ø 4.5		S-BICT45L
Impression Coping Transfer	ø 4.5		S-BICT45S
Impression Coping Transfer	ø 5.5		S-BICT55L
Impression Coping Transfer	ø 5.5		S-BICT55S
Impression Coping Transfer	ø 6.5		S-BICT65L
Impression Coping Transfer	ø 6.5		S-BICT65S
IMPRESSION COPING TRANSFER NON-HEX			
Mini Impression Coping Transfer Non-Hex	ø 3.5		S-BMICTN37L
Mini Impression Coping Transfer Non-Hex	ø 3.5		S-BMICTN37S
Impression Coping Transfer Non-Hex	ø 4.5		S-BICTN45L
Impression Coping Transfer Non-Hex	ø 4.5		S-BICTN45S
Impression Coping Transfer Non-Hex	ø 5.5		S-BICTN55L
Impression Coping Transfer Non-Hex	ø 5.5		S-BICTN55S
Impression Coping Transfer Non-Hex	ø 6.5		S-BICTN65L
Impression Coping Transfer Non-Hex	ø 6.5		S-BICTN65S
FIXTURE LAB ANALOG			
Mini Fixture Lab Analog	ø 3.3	12 mm	S-BMFLA
Fixture Lab Analog	ø 4.3	12 mm	S-BFLA
BONE LEVEL SOLID SYSTEM			
SOLID ABUTMENTS			
Mini Solid Abutment	ø 3.5	1 mm	S-BMSA3710
Mini Solid Abutment	ø 3.5	2 mm	S-BMSA3720
Mini Solid Abutment	ø 3.5	3 mm	S-BMSA3730
Mini Solid Abutment	ø 3.5	4 mm	S-BMSA3740
Mini Solid Abutment	ø 3.5	5 mm	S-BMSA3750

PRODUCT	PLATFORM	HEIGHT	CODE
Solid Abutment	∅ 4.5	1 mm	S-BSA4510
Solid Abutment	∅ 4.5	2 mm	S-BSA4520
Solid Abutment	∅ 4.5	3 mm	S-BSA4530
Solid Abutment	∅ 4.5	4 mm	S-BSA4540
Solid Abutment	∅ 4.5	5 mm	S-BSA4550
Solid Abutment	∅ 5.5	1 mm	S-BSA5510
Solid Abutment	∅ 5.5	2 mm	S-BSA5520
Solid Abutment	∅ 5.5	3 mm	S-BSA5530
Solid Abutment	∅ 5.5	4 mm	S-BSA5540
Solid Abutment	∅ 5.5	5 mm	S-BSA5550
Solid Abutment	∅ 6.5	1 mm	S-BSA6510
Solid Abutment	∅ 6.5	2 mm	S-BSA6520
Solid Abutment	∅ 6.5	3 mm	S-BSA6530
Solid Abutment	∅ 6.5	4 mm	S-BSA6540
Solid Abutment	∅ 6.5	5 mm	S-BSA6550

SOLID PLASTIC BURN-OUT COPING

Mini Solid Plastic Coping Bridge	∅ 3.5		S-BMSAPIC37b
Solid Plastic Coping Bridge	∅ 4.5		S-BSAPIC45b
Solid Plastic Coping Bridge	∅ 5.5		S-BSAPIC55b
Solid Plastic Coping Bridge	∅ 6.5		S-BSAPIC65b
Mini Solid Plastic Coping Crown	∅ 3.5		S-BMSAPIC37c
Solid Plastic Coping Crown	∅ 4.5		S-BSAPIC45c
Solid Plastic Coping Crown	∅ 5.5		S-BSAPIC55c
Solid Plastic Coping Crown	∅ 6.5		S-BSAPIC65c

SOLID ABUTMENT IMPRESSION CAPS

Mini Solid Abutment Impression Cap	∅ 3.5		S-BMSAICap37
Solid Abutment Impression Cap	∅ 4.5		S-BSAICap45
Solid Abutment Impression Cap	∅ 5.5		S-BSAICap55
Solid Abutment Impression Cap	∅ 6.5		S-BSAICap65

SOLID LAB ANALOG

Mini Solid Lab Analog	∅ 3.5		S-BMSALA37
Solid Lab Analog	∅ 4.5		S-BSALA45
Solid Lab Analog	∅ 5.5		S-BSALA55
Solid Lab Analog	∅ 6.5		S-BSALA65

BONE LEVEL OCTA SYSTEM

OCTA ABUTMENTS

Mini Octa Abutment	∅ 3.5	1 mm	S-BMOA10
Mini Octa Abutment	∅ 3.5	2 mm	S-BMOA20
Mini Octa Abutment	∅ 3.5	3 mm	S-BMOA30
Mini Octa Abutment	∅ 3.5	4 mm	S-BMOA40
Mini Octa Abutment	∅ 3.5	5 mm	S-BMOA50

PRODUCT	PLATFORM	HEIGHT	CODE
Octa Abutment	ø 4.8	1 mm	S-BOA10
Octa Abutment	ø 4.8	2 mm	S-BOA20
Octa Abutment	ø 4.8	3 mm	S-BOA30
Octa Abutment	ø 4.8	4 mm	S-BOA40
Octa Abutment	ø 4.8	5 mm	S-BOA50
OCTA ABUTMENT HEALING CAP			
Octa Abutment Healing Cap	ø 4.9	4 mm	S-OAHCap
Octa Abutment Healing Cap Large	ø 4.9	6 mm	S-OAHCAP-L
OCTA SYSTEM PLASTIC CYLINDER OCTA			
Octa System Plastic Cylinder Octa	ø 5.2 mm		S-OAPICy
OCTA SYSTEM PLASTIC CYLINDER NON-OCTA			
Octa System Plastic Cylinder Non-Octa	ø 5.2 mm		S-OAPICyN
OCTA SYSTEM TITANIUM CYLINDER			
Octa System Plastic Cylinder	ø 5.2 mm		S-OATICy
OCTA SYSTEM TITANIUM CYLINDER NON-HEX			
Octa System Plastic Cylinder Non-Hex	ø 5.2		S-OATICyN
OCTA SYSTEM IMPRESSION COPING PICK UP			
Octa System Impression Coping Pick Up	ø 5.5	10 mm	S-ICOAP
OCTA SYSTEM IMPRESSION COPING PICK UP NON-HEX			
Octa System Impression Coping Pick Up Non-Hex	ø 5.5	10 mm	S-ICOAPN
OCTA SYSTEM IMPRESSION COPING TRANSFER			
Octa System Impression Coping Transfer	ø 5.0	7 mm	S-ICOAT
OCTA SYSTEM IMPRESSION COPING TRANSFER NON-HEX			
Octa System Impression Coping Transfer Non-Hex	ø 5.0	7 mm	S-ICOATN
OCTA ABUTMENT LAB ANALOG			
Octa Abutment Lab Analog	ø 4.8	14 mm	S-OALA
BONE LEVEL O-RING SYSTEM			
O-RING ABUTMENTS			
Mini O-Ring Abutment	ø 2.9	0.5 mm	S-BMORA05
Mini O-Ring Abutment	ø 3.5	2 mm	S-BMORA20
Mini O-Ring Abutment	ø 3.5	4 mm	S-BMORA40
O-Ring Abutment	ø 3.42	0.5 mm	S-BORA05
O-Ring Abutment	ø 3.42	1 mm	S-BORA10
O-Ring Abutment	ø 4.5	2 mm	S-BORA20
O-Ring Abutment	ø 4.5	3 mm	S-BORA30
O-Ring Abutment	ø 4.5	4 mm	S-BORA40
O-Ring Abutment	ø 4.5	6 mm	S-BORA60

PRODUCT	PLATFORM	HEIGHT	CODE
RETAINERS			
Retainer	∅ 5.5	2 mm	S-ORR
Retainer	∅ 5.5	3.5 mm	S-ORR-C
Retainer	∅ 5.5	2 mm	S-ORR-O
O-RING LAB ANALOG			
O-Ring Lab Analog	∅ 4.1		S-ORLA
BONE LEVEL LOCATOR SYSTEM			
LOCATOR ABUTMENTS			
Locator Mini Abutment	∅ 3.86mm	1.5mm	S-BMLA15
Locator Mini Abutment	∅ 3.86mm	2.5mm	S-BMLA25
Locator Mini Abutment	∅ 3.86mm	3.5mm	S-BMLA35
Locator Mini Abutment	∅ 3.86mm	4.5mm	S-BMLA45
Locator Mini Abutment	∅ 3.86mm	5.5mm	S-BMLA55
Locator Mini Abutment	∅ 3.86mm	6.5mm	S-BMLA65
Locator Abutment	∅ 3.86mm	1.5mm	S-BLA15
Locator Abutment	∅ 3.86mm	2.5mm	S-BLA25
Locator Abutment	∅ 3.86mm	3.5mm	S-BLA35
Locator Abutment	∅ 3.86mm	4.5mm	S-BLA45
Locator Abutment	∅ 3.86mm	5.5mm	S-BLA55
Locator Abutment	∅ 3.86mm	6.5mm	S-BLA65
LOCATOR MATRIX			
Locator Matrix (includes processing cap)	∅ 5.5mm		S-LMHt
Locator Cap (includes 4 pcs caps)	∅ 4.86 mm		S-LRI
LOCATOR ANALOG			
Locator Lab Analog	∅ 4 mm		S-BLLA
BONE LEVEL SUPERB OVERDENTURE SYSTEM			
SUPERB OVERDENTURE ABUTMENT			
Superb Overdenture Mini Abutment	∅ 2.9 mm	0.5 mm	S-BMEBA05
Superb Overdenture Mini Abutment	∅ 3,5 mm	20 mm	S-BMEBA20
Superb Overdenture Mini Abutment	∅ 3.5mm	40 mm	S-BMEBA40
Superb Overdenture Abutment	∅ 3.42 mm	0.5 mm	S-BEBA05
Superb Overdenture Abutment	∅ 3.42 mm	1 mm	S-BEBA10
Superb Overdenture Abutment	∅ 4.50 mm	2 mm	S-BEBA20
Superb Overdenture Abutment	∅ 4.50 mm	3 mm	S-BEBA30
Superb Overdenture Abutment	∅ 4.50 mm	4 mm	S-BEBA40
Superb Overdenture Abutment	∅ 4.50 mm	5 mm	S-BEBA50
Superb Overdenture Abutment	∅ 4.50 mm	6 mm	S-BEBA60

PRODUCT	PLATFORM	HEIGHT	CODE
SUPERB OVERDENTURE CAP			
Superb Overdenture Cap	ø 5.5 mm	3.5 mm	S-EBOCAP
SUPERB OVERDENTURE LAB ANALOG			
Superb Overdenture Analog	ø 4.1mm	16.6 mm	S-EBLA
BONE LEVEL MULTI-UNIT SYSTEM			
MULTI-UNIT ABUTMENT			
Bone Level Multi Unit Abutment	ø 4.8x17	2.5 mm	S-BMUA48172
Bone Level Multi Unit Abutment	ø 4.8x17	3.5 mm	S-BMUA48173
Bone Level Multi Unit Abutment	ø 4.8x35	2.5 mm	S-BMUA48352
Bone Level Multi Unit Abutment	ø 4.8x35	3.5 mm	S-BMUA48353
MULTI UNIT HEALING CAP			
Multi Unit Healing Cap	ø 4.9	4.0 mm	S-BMUAHCAP
MULTI UNIT IMPRESSION COPING PICK UP			
Bone Level Multi Unit Impression Coping Pick Up	ø 4.8		S-BMUAICP
MULTI UNIT ABUTMENT PLASTIC CYLINDER			
Multi Unit Abutment Plastic Cylinder	ø 4.9		S-MUAPICY
MULTI UNIT LAB ANALOG			
Multi Unit Abutment Lab Analog	ø 4.8	14 mm	S-BMUALA
BAR SYSTEM			
Bar burnout			S-IB-BURNOUT
Bar Positioning			S-IBP-BURNOUT
Bar Housing			S-IBH-BURNOUT
Bar Clip Yellow			S-IBClip-y
Bar Clip Pink			S-IBClip-P
Bar Clip Red			S-IBClip-r
Tool for Inserting Clip			
Key for Parallelometer			
TISSUE LEVEL FIXTURES			
TISSUE LEVEL FIXTURE HYBRID			
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	8 mm	S-TFHR3708
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	10 mm	S-TFHR3710
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	12 mm	S-TFHR3712
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	14 mm	S-TFHR3714
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	8 mm	S-TFHR4308
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	10 mm	S-TFHR4310
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	12 mm	S-TFHR4312
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	14 mm	S-TFHR4314

PRODUCT	PLATFORM	HEIGHT	CODE
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	8 mm	S-TFHR4808
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	10 mm	S-TFHR4810
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	12 mm	S-TFHR4812
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	14 mm	S-TFHR4814
TISSUE LEVEL SHORT FIXTURE STRAIGHT			
Tissue Level Short Fixture Straight SRA	ø 4.3	4 mm	S-TFSR4304
Tissue Level Short Fixture Straight SRA	ø 4.3	5 mm	S-TFSR4305
Tissue Level Short Fixture Straight SRA	ø 4.3	6 mm	S-TFSR4306
Tissue Level Short Fixture Straight SRA	ø 4.8	4 mm	S-TFSR4804
Tissue Level Short Fixture Straight SRA	ø 4.8	5 mm	S-TFSR4805
Tissue Level Short Fixture Straight SRA	ø 4.8	6 mm	S-TFSR4806
Tissue Level Short Fixture Straight SRA	ø 5.5	4 mm	S-TWFSR5504
Tissue Level Short Fixture Straight SRA	ø 5.5	5 mm	S-TWFSR5505
Tissue Level Short Fixture Straight SRA	ø 5.5	6 mm	S-TWFSR5506
TISSUE LEVEL WIDE FIXTURE HYBRID			
Tissue Level Wide Fixture Hybrid SRA	ø 5.5	8 mm	S-BFHR5508
Tissue Level Wide Fixture Hybrid SRA	ø 5.5	10 mm	S-BFHR5510
Tissue Level Wide Fixture Hybrid SRA	ø 5.5	12 mm	S-BFHR5512
TISSUE LEVEL HEALING ABUTMENTS			
Tissue Level Healing Abutment	ø 3.7	3 mm	S-THA30
Tissue Level Healing Abutment	ø 3.7	4 mm	S-THA40
Tissue Level Healing Abutment	ø 3.7	5.5 mm	S-THA55
Tissue Level Healing Abutment	ø 6.4	3 mm	S-THAW30
Tissue Level Healing Abutment	ø 6.4	4 mm	S-THAW40
Tissue Level Healing Abutment	ø 6.4	5.5 mm	S-THAW55
TISSUE LEVEL COUPLE SYSTEM			
TISSUE LEVEL COUPLE ABUTMENTS			
Tissue Level Couple Abutment	ø 5.2	1 mm	S-TCA10
Tissue Level Couple Abutment	ø 5.2	2 mm	S-TCA20
Tissue Level Couple Abutment	ø 5.2	3 mm	S-TCA30
Tissue Level Couple Abutment	ø 5.2	4 mm	S-TCA40
Tissue Level Couple Abutment Wide	ø 6.4	1 mm	S-TCAW10
Tissue Level Couple Abutment Wide	ø 6.4	2 mm	S-TCAW20
Tissue Level Couple Abutment Wide	ø 6.4	3 mm	S-TCAW30
Tissue Level Couple Abutment Wide	ø 6.4	4 mm	S-TCAW40

PRODUCT	PLATFORM	HEIGHT	CODE
TISSUE LEVEL COUPLE ABUTMENT NON-OCTA			
Tissue Level Couple Abutment Non-Octa	ø 5.2	1 mm	S-TCAN10
Tissue Level Couple Abutment Non-Octa	ø 5.2	2 mm	S-TCAN20
Tissue Level Couple Abutment Non-Octa	ø 5.2	3 mm	S-TCAN30
Tissue Level Couple Abutment Non-Octa	ø 5.2	4 mm	S-TCAN40
TISSUE LEVEL COUPLE ABUTMENT WIDE NON-OCTA			
Tissue Level Couple Abutment Non-Octa	ø 6.4	1 mm	S-TCAWN10
Tissue Level Couple Abutment Non-Octa	ø 6.4	2 mm	S-TCAWN20
Tissue Level Couple Abutment Non-Octa	ø 6.4	3 mm	S-TCAWN30
Tissue Level Couple Abutment Non-Octa	ø 6.4	4 mm	S-TCAWN40
TISSUE LEVEL ANGLED ABUTMENTS			
Tissue Level Angled Abutment	ø 3.7	15°	S-TAA15
Tissue Level Angled Abutment	ø 3.7	25°	S-TAA25
TISSUE LEVEL MULTI ABUTMENT OCTA			
Tissue Level Multi Abutment Octa	ø 5.2	1.3 mm	S-TMA
TISSUE LEVEL MULTI ABUTMENT NON-OCTA			
Tissue Level Multi Abutment Non-Octa	ø 5.2	1.3 mm	S-TMAN
TISSUE LEVEL SYNOCTA ABUTMENT OCTA			
Tissue Level Synocta Abutment Octa	ø 3.5	6 mm	S-TSOA
TISSUE LEVEL SYNOCTA ABUTMENT NON-OCTA			
Tissue Level Synocta Abutment Non-Octa	ø 3.5	6 mm	S-TSOAN
TISSUE LEVEL IMPRESSION COPING PICK UP OCTA			
Tissue Level Impression Coping Pick Up Octa	ø 5.5	11.75 mm	S-TICP
Tissue Level Impression Coping Pick Up Wide Octa	ø 6.4	11.75 mm	S-TICPW
TISSUE LEVEL IMPRESSION COPING PICK UP NON-OCTA			
Tissue Level Impression Coping Pick Up Non-Octa	ø 5.5	11.75 mm	S-TICPN
Tissue Level Impression Coping Pick Up Wide Non-Octa	ø 6.4	11.75 mm	S-TICPWN
TISSUE LEVEL IMPRESSION COPING TRANSFER			
Tissue Level Impression Coping Transfer Octa	ø 5.5	9 mm	S-TICT
TISSUE LEVEL IMPRESSION COPING TRANSFER NON-OCTA			
Tissue Level Impression Coping Transfer Non-Octa	ø 5.5	9 mm	S-TICTN
TISSUE LEVEL FIXTURE LAB ANALOG			
Tissue Level Fixture Lab Analog	ø 4.8	14 mm	S-TFLA
Tissue Level Fixture Large Lab Analog	ø 6.0	14 mm	S-TFLLA

PRODUCT	PLATFORM	HEIGHT	CODE
TISSUE LEVEL SOLID SYSTEM			
TISSUE LEVEL SOLID ABUTMENTS			
Tissue Level Solid Abutment	ø 3.5	4 mm	S-TSA40
Tissue Level Solid Abutment	ø 3.5	5.5 mm	S-TSA55
Tissue Level Solid Abutment	ø 3.5	7 mm	S-TSA70
Tissue Level Solid Abutment Large	ø 4.5	4 mm	S-TSAL40
Tissue Level Solid Abutment Large	ø 4.5	5 mm	S-TSAL50
Tissue Level Solid Abutment Large	ø 4.5	7 mm	S-TSAL70
TISSUE LEVEL SOLID ABUTMENT POSITIONING CYLINDER			
Tissue Level Solid Abutment Positioning Cylinder	ø 10.2		S-TSAPoCy
TISSUE LEVEL SOLID ABUTMENT PLASTIC COPING			
Tissue Level Solid Plastic Coping Crown	ø 5.3	10 mm	S-TSAPICc
Tissue Level Solid Plastic Coping Bridge	ø 5.3	10 mm	S-TSAPICb
TISSUE LEVEL SOLID ABUTMENT IMPRESSION CAP			
Tissue Level Solid Abutment Impression Cap	ø 10.2		S-TSAICap
TISSUE LEVEL SOLID ABUTMENT LAB ANALOG			
Tissue Level Solid Abutment Lab Analog	ø 4.8	14 mm	S-TSALA40
Tissue Level Solid Abutment Lab Analog	ø 4.8	15.5 mm	S-TSALA55
Tissue Level Solid Abutment Lab Analog	ø 4.8	17 mm	S-TSALA70
Tissue Level Solid Abutment Large Lab Analog	ø 6.0	14 mm	S-TSALLA40
Tissue Level Solid Abutment Large Lab Analog	ø 6.0	15.5 mm	S-TSALLA55
Tissue Level Solid Abutment Large Lab Analog	ø 6.0	17 mm	S-TSALLA70
TISSUE LEVEL OCTA SYSTEM			
TISSUE LEVEL OCTA ABUTMENT			
Tissue Level Octa Abutment	ø 3.5		S-TOA
OCTA ABUTMENT HEALING CAP			
Octa Abutment Healing Cap	ø 4.9	4 mm	S-OAHCap
Octa Abutment Healing Cap Large	ø 4.9	6 mm	S-OAHCAP-L
TISSUE LEVEL OCTA SYSTEM PLASTIC CYLINDER			
Octa System Plastic Cylinder Octa	ø 5		S-OAPICy
TISSUE LEVEL OCTA ABUTMENT PLASTIC CYLINDER NON-OCTA			
Tissue Level Octa Abutment Plastic Cylinder Non-Octa	ø 5		S-OAPICyN
TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING PICK UP			
Tissue Level Octa System Impression Coping Pick Up Octa	ø 5.5		S-ICOAP
Tissue Level Octa System Impression Coping Pick Up Non-Octa	ø 5.5		S-ICOAPN

PRODUCT	PLATFORM	HEIGHT	CODE
TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING TRANSFER			
Tissue Level Octa System Impression Coping Transfer Octa	∅ 5.0		S-ICOAT
Tissue Level Octa System Impression Coping Transfer Non-Octa	∅ 5.0		S-ICOATN
TISSUE LEVEL OCTA ABUTMENT LAB ANALOG			
Tissue Level Octa Abutment Lab Analog	∅ 4.8		S-OALA
TISSUE LEVEL O-RING SYSTEM			
TISSUE LEVEL O-RING ABUTMENTS			
Tissue Level O-Ring Abutment	∅ 3.5	0.5 mm	S-TORA05
Tissue Level O-Ring Abutment	∅ 3.5	2 mm	S-TORA20
Tissue Level O-Ring Abutment	∅ 3.5	4 mm	S-TORA40
TISSUE LEVEL O-RING LAB ANALOG			
Tissue Level O-Ring Lab Analog	∅ 4.1		S-ORLA
TISSUE LEVEL LOCATOR SYSTEM			
TISSUE LEVEL LOCATOR ABUTMENTS			
Tissue Level Locator Abutment	∅ 4.9	1.5 mm	S-TLA15
Tissue Level Locator Abutment	∅ 4.9	2.5 mm	S-TLA25
Tissue Level Locator Abutment	∅ 4.9	3.5 mm	S-TLA35
Tissue Level Locator Abutment	∅ 4.9	4.5 mm	S-TLA45
Tissue Level Locator Abutment	∅ 4.9	5.5 mm	S-TLA55
Tissue Level Locator Abutment	∅ 4.9	6.5 mm	S-TLA65
LOCATOR MATRIX			
Locator Matrix (includes processing cap)	∅ 5.5mm		S-LMHT
Locator Cap (includes 4 pcs caps)	∅ 4.86 mm		S-LRI
TISSUE LEVEL LOCATOR LAB ANALOG			
Tissue Level Locator Lab Analog			S-TLLA
TISSUE LEVEL SUPERB OVERDENTURE SYSTEM			
TISSUE LEVEL SUPERB OVERDENTURE ABUTMENTS			
Tissue Level Superb Overdenture Abutment	∅ 3.5	0.5 mm	S-TEBA05
Tissue Level Superb Overdenture Abutment	∅ 3.5	2 mm	S-TEBA20
Tissue Level Superb Overdenture Abutment	∅ 3.5	4 mm	S-TEBA40
SUPERB OVERDENTURE CAP			
Superb Overdenture Cap	∅ 5.5 mm	3.5 mm	S-EBOCAP
SUPERB OVERDENTURE LAB ANALOG			
Superb Overdenture Analog	∅ 4.1mm	16.6 mm	S-EBLA

PRODUCT	PLATFORM	HEIGHT	CODE
ZYGOMA IMPLANT			
ABUTMENTS			
Zygoma Abutment	3 mm		S-ZCA4330
Zygoma Abutment	5 mm		S-ZCA4350
ZYGOMA ANGLED ABUTMENTS			
Zygoma Angled Abutment		3 mm	S-ZAA43203
Zygoma Angled Abutment		5 mm	S-ZAA43205
CAPS			
Zygoma Healing Cap		3 mm	S-ZHA4330
Zygoma Healing Cap		5 mm	S-ZHA4350
COPING			
Zygoma Impression Coping Pick Up			S-ZICP
PLASTIC CYLINDER			
Zygoma Plastic Cylinder			S-ZPICYN
LAB ANALOG			
Zygoma Lab Analog			S-ZLA
DRIVERS			
Implant Drivers (For ratchet and contra angle)			
Bone Level Mini Fixture Driver Machine Short			S-BMFDMS
Bone Level Mini Fixture Driver Machine Long			S-BMFDML
Bone Level Mini Fixture Driver Ratchet Short			S-BMFDRS
Bone Level Mini Fixture Driver Ratchet Long			S-BMFDRL
Bone Level Fixture Driver Machine Short			S-BFDMS
Bone Level Fixture Driver Machine Long			S-BFDML
Bone Level Fixture Driver Ratchet Short			S-BFDRS
Bone Level Fixture Driver Ratchet Long			S-BFDRL
Tissue Level Fixture Driver Ratchet Short			S-TFDRS
Tissue Level Fixture Driver Ratchet Long			S-TFDRL
Tissue Level Fixture Driver Machine Short			S-TFDMS
Tissue Level Fixture Driver Machine Long			S-TFDML
Hex driver (extra short, short, long, extra long)			S-HDRS / S-HDRL
Hex Driver Contra Angle (extra short, short, long, extra long)			S-HDRS / S-HDRL

PRODUCT	PLATFORM	HEIGHT	CODE
Superb Overdenture Driver			S-EBRD
O-ring Driver			S-ORADS/L
Octa Driver			S-OADS/L
Solid Driver			S-SADL
ORRO Insert Adaptor			S-OIA
Bar Matrix Insert Adaptor			S-BMIA

SCREWS

CLOSING SCREW

BONE LEVEL

Bone Level Mini Closing Screw	ø 2.8	4.5 mm	S-BMCS
Bone Level Closing Screw	ø 3.36	6.5 mm	S-BCS

TISSUE LEVEL

Tissue Level Closing-Screw Small	ø 3.5	5.35 mm	S-TCSS
Tissue Level Closing-Screw Large	ø 5.2	6.4mm	S-TCSL
Tissue Level Closing-Screw Wide	ø 6.4	6.4 mm	S-TCSW

ABUTMENT SCREWS

BONE LEVEL

Bone Level Mini Abutment Screw	ø 2.1	7.64 mm	S-BMAS
Bone Level Abutment Screw	ø 2.3	9.32 mm	S-BAS
Bone Level Multi-Unit Abutment Screw	ø 1.96	9 mm	S-BMUAS

TISSUE LEVEL

Tissue Level Abutment Screw	ø 2.5	7.1	S-TAS
-----------------------------	-------	-----	-------

COMMON SCREWS

Octa Abutment Screw	ø 2.5	4.5 mm	S-OAS
Multi-Unit Abutment Healing Screw	ø 2.2	4 mm	S-MUAHS

Blank lined writing area consisting of 30 horizontal lines for notes.

